

Facts:

- One out of three Americans over the age of 65 fall at home each year.
- Most home accidents are compounded because the injured person lays unattended for an extended period.
- The chances of surviving a fall are six times greater if you are attended to within one hour.
- Being unable to get help when needed is a common event that usually marks the end of a person's independent living.

HouseCall

Personal Emergency Response System




A Key to Independent Living

We at **Lighthouse Systems** have over 28 years of alarm experience to offer our customers.

Our products focus on one very fundamental premise: **providing our customers the highest level of peace of mind by combining the best products, services and support available.**

We'll show you the way...



...to greater peace of mind.

Lighthouse Systems, LLC

3820 S Jackson Rd, Unit A
Pharr, TX 78577
Phone: 956-631-7016
Toll Free: 1-866-232-0415
Fax: 956-631-7016
E-mail: info@lighthouse-systems.com

Texas Department of State Health Services

PERS License # PP175

Texas Department of Aging and Disability Services


CCADIERS Vendor # 1014703

CBAIERS Vendor # 1014702

Lighthouse Systems, LLC

HouseCall

Personal Emergency Response System



**A Key to Peace of Mind
and Independent Living**



1-866-232-0415



Lighthouse Systems, LLC

"We'll show you the way...
...To greater peace of mind."

Lighthouse Systems, LLC

HouseCall

Personal Emergency Response System

Help when you need it most!

We can't always prevent an emergency from occurring. Because household accidents do happen, it's important that



Help is just a push of a button away.

we get help fast when needed. In many cases, it's the delay in getting help that creates a debilitating problem.

With a simple push of a button,

HouseCall can have the help you need on the way. You don't even have to get

to the phone! A trained dispatcher will talk to you via the main medical alarm console. They will stay with you until the help you need arrives. You determine the



Transmitters available in a watch or pendent

assistance you receive. In case you are unable to respond to the dispatcher, we will follow the instructions you have previously provided.

HouseCall Features:

- **Two-way voice communication:** talk back and forth with our trained emergency dispatcher.
- **Large emergency push button on main console:** it's easy to find in an emergency.
- **Portable emergency transmitter:** they come in pendant or watch styles, are water-resistant (can be worn in the bath) and transmit over 100 feet.
- **Braille characters for the main function buttons on console:** it's easy to use for the sight impaired.
- **Selectable inactivity reporting: 8, 12 or 24 hours (optional):** if there is no activity in the house, we will call and check up on you.
- **Speakerphone ability:** push the button on your emergency transmitter to have the console answer your incoming phone call.
- **Capacity for up to 16 emergency transmitters:** each person can have their own unit and have additional transmitters placed in strategic, high risk areas, such as the bathroom.
- **32-hour battery backup:** the system is still there for you even if the power is out.

"...Hello, This is the emergency monitoring center. Is everything okay?..."




Local Company, Nation-wide Service

HouseCall

Personal Emergency Response System
Take it with you when you relocate or travel.

The Lighthouse Medical Alert System can operate wherever there is electrical power and an active telephone line. The unit is lightweight and completely portable. **Re-installing and testing the unit takes less than five minutes!**

To update your contact information, simply call the toll free number listed on the top of the main console. It's that easy!



Call us now for more information:

956-631-7016

Toll Free

1-866-232-0415

We'll show you the way...



Lighthouse Systems, LLC

3820 S Jackson Rd, Unit A
Pharr, TX 78577
Phone: 956-631-7016
Toll Free: 1-866-232-0415
Fax: 956-631-7016

E-mail: info@lighthouse-systems.com

Texas Department of State Health Services
PERS License # PPI75

Texas Department of Aging and Disability Services
CCAD/ERS Vendor # 1014703
CBA/ERS Vendor # 1014702

www.lighthouse-systems.com